

PLDAPS syncing

- Multiple clocks/timings need to be synchronized
 - PTB computer :: Datapixx box
 - PDS data :: Ephys (e.g. Plexon OmniPlex/MAP, OpenEphys, ...)
 - PTB fliptime :: Stimulus display onset
 - PTB computer :: Eyelink
- Test & calibration code should explain where in your rigPrefs the results & settings for these syncs should be stored so that they can be applied/paired with your data files by default
 - e.g. `createRigPrefs(additionalSettings)`

PTB :: Datapixx

`PsychDataPixx('GetPreciseTime')`

- Prevents drift between PTB (GetSecs.m) & Datapixx clocks
- Computational cost (≥ 20 ms) too large to be done during active trials
- Called automatically at the start of every trial
by: `pldapsDefaultTrialFunction.m >> trialSetup`

`PLDAPS/testsuite/pldapsDatapixxSync.m`

- Test code for selecting optimal low-level parameters for your system
- Describes how/where to save those values in your RigPrefs